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Ex Classification  Ex ia IIC T6/T5/T4/T3 Ga/Gb in accordance to the test certificate number SEV 14 ATEX 0168 X; IECEx SEV 14.0025X and CML 22 UKEX 2108X

1 Introduction

According to Directive 2014/34/EU (ATEX 114)¹⁾ Appendix I, the pH sensors type InPro™ 2XXX, InPro 3XXX, InPro 4XXX and InPro X1 are devices of equipment group II; category 1/2G which, according to Directive 99/92/EC (ATEX 137)²⁾ can be used in zones 0/1 or 1/2 or 1 or 2 as well as gas groups IIA, IIB and IIC which are potentially explosive due to combustible substances in temperature classes T3 to T6.

The requirements specified in EN 60079-14 must be observed during use/installation.

For the analog version of the pH electrode, the pH measurement circuit, temperature measurement circuit and data chip circuit are part of a common intrinsically safe system and are jointly connected to and operated by a separately certified transmitter.

The digital version of the pH sensor is connected to and operated by a two-wire cable to the certified transmitter.

The intrinsically safe circuits are galvanically isolated from the non-intrinsically safe circuits up to a nominal voltage peak value of 375 V and from the earthed parts up to a nominal voltage peak value of 30 V.

2 Specifications

Analog pH sensor	With type of protection: intrinsic safety to Ex ia IIC
pH measuring circuit, temperature measuring circuit and data chip circuit	Only for connection to certified intrinsically safe circuits.
Maximum values	$U_i \leq 16\text{V}$, $I_i \leq 30\text{ mA}$, $P_i \leq 50\text{ mW}$ or $U_i \leq 16\text{V}$, $I_i \leq 190\text{ mA}$, $P_i \leq 200\text{ mW}$ $L_i = 0$ (effective internal inductance) $C_i = 900\text{ pF}$ (effective internal capacitance)

In each case, the above values shall represent the sum of the individual circuits in the intrinsically safe power supply and evaluation device concerned (transmitter).

Digital pH sensor	With type of protection: intrinsic safety to Ex ia IIC
Two-wire current circuit	Only for connection to certified intrinsically safe circuits.
Maximum values	$U_i \leq 16\text{V}$, $I_i \leq 30\text{ mA}$, $P_i \leq 50\text{ mW}$ $L_i = \text{negligible}$ $C_i = \text{negligible}$

¹⁾ For UK statutory Requirements SI 2016 No. 1107

²⁾ For UK regulations 7 and 11 of the Dangerous Substances and Explosive Atmospheres Regulations 2002 (DSEAR)

3 Special conditions

3.1 The relationship between the maximum permissible ambient or media temperature and temperature class is shown in the following table:

With analog pH sensor: For $U_i \leq 16\text{ V}$, $I_i \leq 30\text{ mA}$, $P_i \leq 50\text{ mW}$; pH measuring circuit, temperature measurement circuit, and data chip circuit:	Temperature class	Maximum ambient or media temperature
	T6	62 °C
	T5	74 °C
	T4	102 °C
	T3	154 °C

or

With analog pH sensor: For $U_i \leq 16\text{ V}$, $I_i \leq 190\text{ mA}$, $P_i \leq 200\text{ mW}$; pH measuring circuit, temperature measurement circuit, and data chip circuit:	Temperature class	Maximum ambient or media temperature
	T6	51 °C
	T5	63 °C
	T4	91 °C
	T3	143 °C

or

With digital pH sensor: For $U_i \leq 16\text{ V}$, $I_i \leq 30\text{ mA}$, $P_i \leq 50\text{ mW}$; two-wire current circuit:	Temperature class	Maximum ambient or media temperature
	T6	62 °C
	T5	74 °C
	T4	102 °C
	T3	131 °C

3.2 The pH sensor types InPro 327X, InPro 427X, InPro 487X and InPro X1 are constructed from plastic. To prevent the risk of electrostatic sparking, the plastic surface should only be cleaned with a damp cloth.

3.3 The enclosures of pH sensor types InPro 328X, InPro 428X, InPro 488X and InPro X1 containing titanium constitute a potential risk of ignition by impact or friction. Care must be taken into account during installation and use to prevent impact or friction.

3.4 For safe working and to prevent explosion, the minimum conductivity of media must be higher than 1 nS/cm when the pH sensor types of InPro 327X, InPro 427X, InPro 487X and InPro X1 are used.

3.5 The capacitance and inductance of the connecting cable has to be considered.

3.6 The pH sensors types InPro 2XXX, InPro 3XXX, InPro 4XXX and InPro X1 can be used in/with the fittings InFit™ 76X-XXX or InTrac™ 7XX-XXX, or in/with other suitable fittings in potentially explosive areas.

The metal body of the pH sensors, or the fitting InFit 76X-XXX or InTrac 7XX-XXX, or other appropriate fitting is optionally included in the routine pressure test of the system.

3.7 The independent fitting used for installation of pH electrodes must be conductively connected to the equipotential bonding system.

InPro, InFit and InTrac are trademarks of the METTLER TOLEDO Group.

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All rights reserved